

KIF16B

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:152kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human KIF16B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SNX23; C20orf23; KISC20ORF;

Background:

The protein encoded by this gene is a kinesin-like protein that may be involved in intracellular trafficking. Three transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).

Catalog #:A10314

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:55614

Isotype:IgG

Swiss Prot:Q96L93

Purity:Affinity purification

For research use only.