

## CISD1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:1000 - 1:2000

**Calculated MW:**12kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human CISD1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ZCD1; MDS029; C10orf70; mitoNEET;

**Background:**

This gene encodes a protein with a CDGSH iron-sulfur domain and has been shown to bind a redox-active

*To place an order, please [Click HERE](#).*

**Catalog #:**A10317

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**55847

**Isotype:**IgG

**Swiss Prot:**Q9NZ45

**Purity:**Affinity purification

For research use only.