

## GPR126

---

**Reactivity:**Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**137kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human GPR126

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

APG1; DREG; VIGR; LCCS9; GPR126; PS1TP2;

**Background:**

This gene, which is upregulated in human umbilical vein endothelial cells, encodes a G protein-coupled receptor. Variations in this gene can affect a person's stature. Multiple transcript variants encoding different proteins have been found for this gene.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10318

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**57211

**Isotype:**IgG

**Swiss Prot:**Q86SQ4

**Purity:**Affinity purification

For research use only.