

LRRC4

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:73kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human LRRC4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NAG14; NGL-2;

Background:

This gene is significantly downregulated in primary brain tumors. The exact function of the protein encoded by this gene is unknown.

To place an order, please [Click HERE](#).

Catalog #:A10321

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:64101

Isotype:IgG

Swiss Prot:Q9HBW1

Purity:Affinity purification

For research use only.