

SPG11

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:279kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SPG11

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ALS5; CMT2X; KIAA1840;

Catalog #:A10326

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:80208

Isotype:IgG

Swiss Prot:Q96J17

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a potential transmembrane protein that is phosphorylated upon DNA damage. Defects in this gene are a cause of spastic paraplegia type 11 (SPG11). Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).