

JPH4

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 66kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human JPH4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

JP4; JPHL1;

Catalog #: A10333

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84502

Isotype: IgG

Swiss Prot: Q96JJ6

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the junctophilin family of transmembrane proteins that are involved in the formation of the junctional membrane complexes between the plasma membrane and the endoplasmic/sarcoplasmic reticulum in excitable cells. The encoded protein contains a conserved N-terminal repeat region called the membrane occupation and recognition nexus sequence that is found in other members of the junctophilin family. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).