

STRIP1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:96kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human STRIP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FAM40A; FAR11A;

Catalog #:A10334

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:85369

Isotype:IgG

Swiss Prot:Q5VSL9

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the striatin-interacting phosphatase and kinase complex, which is involved in localization of the Golgi body. The encoded protein participates in cytoskeletal organization. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).