

EGLN1/EGLN2

Reactivity:Human

Tested applications:WB

Calculated MW:426aa/407aa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human EGLN1/EGLN2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10342

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:54583/112398

Isotype:IgG

Swiss Prot:Q9GZT9/Q96KS0

Purity:Affinity purification

For research use only.