

## CYP2R1

---

**Reactivity:**Mouse Rat

**Tested applications:**WB

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human CYP2R1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**To place an order, please [Click HERE](#).**

**Catalog #:**A10470

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**120227

**Isotype:**IgG

**Swiss Prot:**Q6VVX0

**Purity:**Affinity purification

For research use only.