

## PRLHR

---

**Reactivity:**Human Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:1000

**Calculated MW:**41kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human PRLHR

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

GR3; GPR10; PrRPR;

**Background:**

PRLHR is a 7-transmembrane domain receptor for prolactin-releasing hormone (PRLH; MIM 602663) that is highly expressed in anterior pituitary (Ozawa et al., 2002

*To place an order, please [Click HERE](#).*

**Catalog #:**A10479

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**2834

**Isotype:**IgG

**Swiss Prot:**P49683

**Purity:**Affinity purification

For research use only.