

AKR1C2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:37kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human AKR1C2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AKR1C2;AKR1C-pseudo;BABP;DD;DD2;DDH2;HAKRD;HBAB;MCDR2 ;

Background:

This gene encodes a member of the aldo/keto reductase superfamily, which consists of more than 40 known enzymes and proteins. These enzymes catalyze the conversion of aldehydes and ketones to their corresponding alcohols using NADH and/or NADPH as cofactors. The enzymes display overlapping but distinct substrate specificity. This enzyme binds bile acid with high affinity, and shows minimal 3-alpha-hydroxysteroid dehydrogenase activity. This gene shares high sequence identity with three other gene members and is clustered with those three genes at chromosome 10p15-p14. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Dec 2011]

To place an order, please [Click HERE](#).

Catalog #:A1048

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1646

Isotype:IgG

Swiss Prot:P52895

Purity:Affinity purification

For research use only.