

SLC22A4

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:62kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SLC22A4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OCTN1;

Background:

Polyspecific organic cation transporters in the liver, kidney, intestine, and other organs are critical for elimination of many endogenous small organic cations as well as a wide array of drugs and environmental toxins. The encoded protein is an organic cation transporter and plasma integral membrane protein containing eleven putative transmembrane domains as well as a nucleotide-binding site motif. Transport by this protein is at least partially ATP-dependent.

To place an order, please [Click HERE](#).

Catalog #:A10490

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6583

Isotype:IgG

Swiss Prot:Q9H015

Purity:Affinity purification

For research use only.