

SLC22A4

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 62kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SLC22A4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OCTN1;

Catalog #: A10490

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6583

Isotype: IgG

Swiss Prot: Q9H015

Purity: Affinity purification

For research use only.

Background:

Polyspecific organic cation transporters in the liver, kidney, intestine, and other organs are critical for elimination of many endogenous small organic cations as well as a wide array of drugs and environmental toxins. The encoded protein is an organic cation transporter and plasma integral membrane protein containing eleven putative transmembrane domains as well as a nucleotide-binding site motif. Transport by this protein is at least partially ATP-dependent.

To place an order, please [Click HERE](#).