

DSG4

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:114kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human DSG4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LAH; HYPT6; CDGF13; CDHF13

Catalog #:A10502

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:147409

Isotype:IgG

Swiss Prot:Q86SJ6

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the desmoglein subgroup of desmosomal cadherins. The encoded protein is a transmembrane component in desmosomes and may play a role in cell-cell adhesion in epithelial cells. Mutations in the gene are associated with localized autosomal recessive hypotrichosis and potentially in other skin disorders. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).