

PCDH17

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:126kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PCDH17

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PCH68; PCDH68;

Background:

This gene belongs to the protocadherin gene family, a subfamily of the cadherin superfamily. The encoded protein contains six extracellular cadherin domains, a transmembrane domain, and a cytoplasmic tail differing from those of the classical cadherins. The encoded protein may play a role in the establishment and function of specific cell-cell connections in the brain.

To place an order, please [Click HERE](#).

Catalog #:A10512

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:27253

Isotype:IgG

Swiss Prot:O14917

Purity:Affinity purification

For research use only.