

## C16orf54

---

**Reactivity:**Mouse Rat

**Tested applications:**WB

**Calculated MW:**225 aa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human C16orf54

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10519

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**105369277/283897

**Isotype:**IgG

**Swiss Prot:**Q8C708

**Purity:**Affinity purification

For research use only.