

Hisne H3 (tri methyl K27)

Reactivity:Human

Tested applications:WB IF

Observed MW:Refer to figures

Immunogen:

Synthetic Peptide

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10529

Antibody Type:

Monoclonal Antibody

Species:Mouse

Gene ID:8290

Isotype:IgG

Swiss Prot:Q16695

Purity:Affinity purification

For research use only.