

UBE2N

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 17kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human UBE2N

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UBC13; Ubch13; HEL-S-71; Ubch-ben; UBCHBEN; UBC13;

Catalog #: A10543

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7334

Isotype: IgG

Swiss Prot: P61088

Purity: Affinity purification

For research use only.

Background:

The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. Studies in mouse suggest that this protein plays a role in DNA postreplication repair.

To place an order, please [Click HERE](#).