

## GIPC1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:1000 - 1:2000

**Calculated MW:**36kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human GIPC1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

NIP; GIPC; IIP-1; TIP-2; SEMCAP; C19orf3; Hs.6454; GLUT1CBP; RGS19IP1; SYNECTIN;  
SYNECTIIN;

**Background:**

GIPC1 is a scaffolding protein that regulates cell surface receptor expression and trafficking (Lee et al., 2008

*To place an order, please [Click HERE](#).*

**Catalog #:**A10550

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**10755

**Isotype:**IgG

**Swiss Prot:**O14908

**Purity:**Affinity purification

For research use only.