

SEC11A

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:21kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SEC11A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SPC18; SPCS4A; SEC11L1; sid2895; 1810012E07Rik;

Catalog #:A10552

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23478

Isotype:IgG

Swiss Prot:P67812

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the peptidase S26B family. The encoded protein is an 18kDa subunit of the signal peptidase complex and has been linked to cell migration and invasion, gastric cancer and lymph node metastasis. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 8.

To place an order, please [Click HERE](#).