

AP1S1

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:19kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human AP1S1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AP19; EKV3; CLAPS1; MEDNIK; SIGMA1A;

Catalog #:A10558

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1174

Isotype:IgG

Swiss Prot:P61966

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is part of the clathrin coat assembly complex which links clathrin to receptors in coated vesicles. These vesicles are involved in endocytosis and Golgi processing.

This protein, as well as beta-prime-adaptin, gamma-adaptin, and the medium (mu) chain AP47, form the AP-1 assembly protein complex located at the Golgi vesicle.

To place an order, please [Click HERE](#).