

PARG

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:111kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PARG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PARG99;

Catalog #:A10577

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8505

Isotype:IgG

Swiss Prot:Q86W56

Purity:Affinity purification

For research use only.

Background:

Poly(ADP-ribose) glycohydrolase (PARG) is the major enzyme responsible for the catabolism of poly(ADP-ribose), a reversible covalent-modifier of chromosomal proteins. The protein is found in many tissues and may be subject to proteolysis generating smaller, active products. Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).