

SNF8

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:29kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SNF8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Dot3; EAP30; VPS22;

Catalog #:A10583

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:11267

Isotype:IgG

Swiss Prot:Q96H20

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a component of the endosomal sorting complex required for transport II (ESCRT-II), which regulates the movement of ubiquitinated transmembrane proteins to the lysosome for degradation. This complex also interacts with the RNA polymerase II elongation factor (ELL) to overcome the repressive effects of ELL on RNA polymerase II activity. Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).