

USP24

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:294kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human USP24

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A10585

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23358

Isotype:IgG

Swiss Prot:Q9UPU5

Purity:Affinity purification

For research use only.

Background:

Modification of cellular proteins by ubiquitin is an essential regulatory mechanism controlled by the coordinated action of multiple ubiquitin-conjugating and deubiquitinating enzymes. USP24 belongs to a large family of cysteine proteases that function as deubiquitinating enzymes (Quesada et al., 2004

To place an order, please [Click HERE](#).