

WDR26

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 72kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human WDR26

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CDW2; GID7; MIP2;

Catalog #: A10594

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 80232

Isotype: IgG

Swiss Prot: Q9H7D7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-aspartate (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. Two transcript variants encoding two different isoforms have been found for this gene.

To place an order, please [Click HERE](#).