

TRMT61A

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 31kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TRMT61A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GCD14; TRM61; Gcd14p; hTRM61; C14orf172;

To place an order, please [Click HERE](#).

Catalog #: A10597

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 115708

Isotype: IgG

Swiss Prot: Q96FX7

Purity: Affinity purification

For research use only.