

WDFY2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 45kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human WDFY2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PROF; WDF2; ZFYVE22;

Catalog #: A10598

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 115825

Isotype: IgG

Swiss Prot: Q96P53

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that contains two WD domains and an FYVE zinc finger region. The function of this gene is unknown. An alternatively spliced transcript variant of this gene may exist.

To place an order, please [Click HERE](#).