

## VKORC1L1

---

**Reactivity:**Rat

**Tested applications:**WB

**Calculated MW:**176aa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human VKORC1L1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10600

**Antibody Type:**

Polyclonal Antibody

**Species:**Rat

**Gene ID:**154807

**Isotype:**IgG

**Swiss Prot:**Q6TEK5

**Purity:**Affinity purification

For research use only.