

## GPHA2

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:1000

**Calculated MW:**14kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human GPHA2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

A2; GPA2; ZSIG51;

**Background:**

GPHA2 is a cystine knot-forming polypeptide and a subunit of the dimeric glycoprotein hormone family (Hsu et al., 2002)

*To place an order, please [Click HERE](#).*

**Catalog #:**A10601

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**170589

**Isotype:**IgG

**Swiss Prot:**Q96T91

**Purity:**Affinity purification

For research use only.