

LYZL6

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:17kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human LYZL6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LYZL6;1700023H08Rik;LYC1;PRO1485;TKAL754 ;

Catalog #:A1067

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:57151

Isotype:IgG

Swiss Prot:O75951

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the C-type lysozyme/alpha-lactalbumin family. C-type lysozymes are bacteriolytic factors that play a role in host defense, whereas alpha-lactalbumins mediate lactose biosynthesis. The encoded protein contains catalytic residues characteristic of C-type lysozymes and may play a role in male reproduction. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jan 2011]

To place an order, please [Click HERE](#).