

## AKR1A1

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

**Calculated MW:**37kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human AKR1A1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

AKR1A1;ALDR1;ALR;ARM;DD3;MGC12529;MGC1380 ;

**Background:**

This gene encodes a member of the aldo/keto reductase superfamily, which consists of more than 40 known enzymes and proteins. This member, also known as aldehyde reductase, is involved in the reduction of biogenic and xenobiotic aldehydes and is present in virtually every tissue. Multiple alternatively spliced transcript variants of this gene exist, all encoding the same protein. [provided by RefSeq, Jan 2011]

*To place an order, please [Click HERE](#).*

**Catalog #:**A1068

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**10327

**Isotype:**IgG

**Swiss Prot:**P14550

**Purity:**Affinity purification

For research use only.