

## KCNJ2

---

**Reactivity:**Human Mouse

**Tested applications:**WB

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human KCNJ2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**To place an order, please [Click HERE](#).**

**Catalog #:**A10708

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**3759

**Isotype:**IgG

**Swiss Prot:**P63252

**Purity:**Affinity purification

For research use only.