

AQP1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50- 1:100

Calculated MW:29kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human AQP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CO; CHIP28; AQP-CHIP

Catalog #:A10730

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:358

Isotype:IgG

Swiss Prot:P29972

Purity:Affinity purification

For research use only.

Background:

Aquaporins are a family of small integral membrane proteins related to the major intrinsic protein (MIP or AQP0). This gene encodes an aquaporin which functions as a molecular water channel protein. It is a homotetramer with 6 bilayer spanning domains and N-glycosylation sites. The protein physically resembles channel proteins and is abundant in erythrocytes and renal tubes. The gene encoding this aquaporin is a possible candidate for disorders involving imbalance in ocular fluid movement. Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).