

## DPYSL5

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Calculated MW:**564

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human DPYSL5

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10734

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**56896

**Isotype:**IgG

**Swiss Prot:**Q9BPU6

**Purity:**Affinity purification

For research use only.