

## POLR3A

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**156kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human POLR3A

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ADDH; HLD7; RPC1; RPC155; hRPC155

**Catalog #:**A10737

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**11128

**Isotype:**IgG

**Swiss Prot:**O14802

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is the catalytic component of RNA polymerase III, which synthesizes small RNAs. The encoded protein also acts as a sensor to detect foreign DNA and trigger an innate immune response.

*To place an order, please [Click HERE](#).*