

IL7R

Reactivity:Human Mouse Rat

Tested applications:WB IHC ICC

Calculated MW:459

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human IL7R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10740

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3575

Isotype:IgG

Swiss Prot:P16871

Purity:Affinity purification

For research use only.