

## HSP90B1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC ICC

**Calculated MW:**803

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human HSP90B1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10746

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**7184

**Isotype:**IgG

**Swiss Prot:**P14625

**Purity:**Affinity purification

For research use only.