

DiMethyl-Histone H3-K4

Reactivity:Human Mouse Rat

Tested applications:WB IHC ICC FC

Calculated MW:136

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human DiMethyl-Histone H3-K4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10748

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8350

Isotype:IgG

Swiss Prot:P68431

Purity:Affinity purification

For research use only.