

HAPLN1

Reactivity:Human Mouse Rat

Tested applications:WB IHC ICC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:48 kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human HAPLN1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CRT1; CRTL1;

To place an order, please [Click HERE](#).

Catalog #:A10828

Antibody Type:

Monoclonal Antibody

Species:Rabbit

Gene ID:1404

Isotype:IgG

Swiss Prot:P10915

Purity:Affinity purification

For research use only.