

## FOLH1

---

**Reactivity:**Human

**Tested applications:**WB IHC IF

**Calculated MW:**750

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human FOLH1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10863

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**2346

**Isotype:**IgG

**Swiss Prot:**Q04609

**Purity:**Affinity purification

For research use only.