

## IBSP

---

**Reactivity:**Human

**Tested applications:**WB

**Calculated MW:**317

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human IBSP

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**To place an order, please [Click HERE](#).**

**Catalog #:**A10867

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**3381

**Isotype:**IgG

**Swiss Prot:**P21815

**Purity:**Affinity purification

For research use only.