

## SRP72

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF FC

**Recommended Dilution:**WB 1:500 - 1:1000 IHC 1:50- 1:100 IF 1:20- 1:50 FC 1:20- 1:50

**Calculated MW:**75kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human SRP72

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10899

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**6731

**Isotype:**IgG

**Swiss Prot:**O76094

**Purity:**Affinity purification

For research use only.