

ICAM1

Reactivity:Human

Tested applications:WB IHC ICC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:89 kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ICAM1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BB2; CD54; P3.58;

Background:

This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor.

To place an order, please [Click HERE](#).

Catalog #:A10907

Antibody Type:

Monoclonal Antibody

Species:Rabbit

Gene ID:3383

Isotype:IgG

Swiss Prot:P05362

Purity:Affinity purification

For research use only.