

## IFNAR1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IP FC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IP 1:20 - 1:50 FC 1:20 - 1:50

**Calculated MW:**90/130 kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human IFNAR1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

AVP; IFRC; IFNAR; IFNBR; IFN-alpha-REC;

**Background:**

The protein encoded by this gene is a type I membrane protein that forms one of the two chains of a receptor for interferons alpha and beta. Binding and activation of the receptor stimulates Janus protein kinases, which in turn phosphorylate several proteins, including STAT1 and STAT2. The encoded protein also functions as an antiviral factor.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10913

**Antibody Type:**

Monoclonal Antibody

**Species:**Rabbit

**Gene ID:**3454

**Isotype:**IgG

**Swiss Prot:**P17181

**Purity:**Affinity purification

For research use only.