

ARID5B

Reactivity:Human Rat

Tested applications:WB

Calculated MW:1188

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human ARID5B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10950

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:84159

Isotype:IgG

Swiss Prot:Q14865

Purity:Affinity purification

For research use only.