

Phospho-PTK2-S722

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF IP

Observed MW:Refer to figures

Immunogen:

Synthetic Peptide

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10958

Antibody Type:

Monoclonal Antibody

Species:Mouse

Gene ID:5747

Isotype:IgG

Swiss Prot:Q05397

Purity:Affinity purification

For research use only.