

Phospho-CAMK2A-T286

Reactivity:Rat

Tested applications:WB IF IP

Observed MW:Refer to figures

Immunogen:

Synthetic Peptide

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10970

Antibody Type:

Monoclonal Antibody

Species:Mouse

Gene ID:815

Isotype:IgG

Swiss Prot:Q9UQM7

Purity:Affinity purification

For research use only.