

## Phospho-p44/42 MAPK -T202/Y204

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IF IP

**Observed MW:**Refer to figures

**Immunogen:**

Synthetic Peptide

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10976

**Antibody Type:**

Monoclonal Antibody

**Species:**Mouse

**Gene ID:**5594/5595

**Isotype:**IgG

**Swiss Prot:**P28482

**Purity:**Affinity purification

For research use only.