

## FGL1

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50- 1:200

**Calculated MW:**36kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic Peptide of human FGL1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HFREP1; HP-041; LFIRE1; LFIRE-1

**Catalog #:**A10977

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**2267

**Isotype:**IgG

**Swiss Prot:**Q08830

**Purity:**Affinity purification

For research use only.

**Background:**

Fibrinogen-like 1 is a member of the fibrinogen family. This protein is homologous to the carboxy terminus of the fibrinogen beta- and gamma- subunits which contains the four conserved cysteines of fibrinogens and fibrinogen related proteins. However, this protein lacks the platelet-binding site, cross-linking region and a thrombin-sensitive site which are necessary for fibrin clot formation. This protein may play a role in the development of hepatocellular carcinomas. Four alternatively spliced transcript variants encoding the same protein exist for this gene.

*To place an order, please [Click HERE](#).*