

ATP1A1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF IP

Calculated MW:1023

Observed MW:Refer to figures

Immunogen:

Recombinant protein ()

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10983

Antibody Type:

Monoclonal Antibody

Species:Mouse

Gene ID:476

Isotype:IgG

Swiss Prot:P05023

Purity:Affinity purification

For research use only.